


Mag-Bind® Viral RNA Xpress Kit

High throughput isolation of viral RNA from nasopharyngeal swab specimens



Isolation from swabs

Reliable isolation of viral RNA from NP swabs that are dry or in VTM



No Pro K Required

Process faster with no Proteinase K digestion step



Improved RT-qPCR Performance

Efficient recovery of viral nucleic acids with minimal inhibitor carryover for demanding applications such as RT-qPCR



Bead-Based

High throughput, scalable purification



Automatable

Adaptable to most liquid handling platforms

Mag-Bind® Viral RNA Express Kit follows a magnetic bead-based approach for the rapid and reliable isolation of viral RNA from nasopharyngeal (NP) swab specimens that are dry or in viral transport media (VTM). The extraction methodology is easily adaptable to various sample types, automated systems and can also be scaled up or down depending on the amount of starting sample amount used. The kit utilizes the proven Mag-Bind® technology that enables purification of high-quality nucleic acids that are free of proteins, nucleases, and other impurities. The purified nucleic acids are ready for direct use in downstream applications such as qPCR, RT-qPCR and more.

Detection of Synthetic SARS-CoV-2 virus control following RT-qPCR




Figure 1. 1×10^5 copies of synthetic SARS-CoV-2 was spiked into a 200 µL sample containing 2000 HEK293 cells. Viral nucleic acids were extracted following the recommended protocol from Mag-Bind® Viral RNA Xpress Kit (M6219). 4 and 8 µL of template was used in a 20 µL SYBR Green-labeled RT-qPCR reaction. The average Ct values obtained are shown on the left. The Ct difference between the two template amounts is ~1 indicating no qPCR inhibition.

Detection of Influenza B virus following RT-qPCR: Mag-Bind® Viral RNA Xpress Kit vs Mag-Bind® Viral DNA/RNA 96 Kit




Figure 2. 50 µL of Influenza B virus control was spiked into a 200 µL sample containing 2000 HEK293 cells. Viral nucleic acids were extracted following the recommended protocols from Mag-Bind® Viral RNA Xpress Kit (M6219) and Mag-Bind® Viral DNA/RNA 96 Kit (M6246). 4 and 8 µL of template was used in a 20 µL SYBR Green-labeled RT-qPCR reaction. The average Ct values obtained are shown left. The Ct difference between the two template amounts used is ~1 indicating no qPCR inhibition. The Ct values for the Mag-Bind® Viral RNA Xpress Kit (M6219) were on an average ~0.5 Ct lower than Mag-Bind® Viral DNA/RNA 96 Kit (M6246) indicating improved performance.

Automation protocols available for:

- Microlab® STAR
- Microlab® MagEx STARlet
- Microlab® NIMBUS
- KingFisher™, BioSprint®, MagMax® 96
- Tecan DreamPrep™ NAP

Product Description	Cat No.
Mag-Bind® Viral RNA Xpress Kit 24x96 preps	M6219-2304



innovations in nucleic acid isolation

Omega Bio-tek, Inc.
400 Pinnacle Way, Suite 450
Norcross, GA 30071

Phone: 770-931-8400
Email: info@omegabiotek.com
Web: www.omegabiotek.com



© 2020 Omega Bio-tek, Inc. All rights reserved. Mag-Bind is a trademark of Omega Bio-tek, Inc. All other trademarks and trade names are the property of their respective holders. For research use only.

Lit No. SL-0125